



Advanced Pharmacy Fridges



RPFG44043 Glass Door Advanced Pharmacy Refrigerator



CLAIM
YOUR
FREE **2**
YEAR
WARRANTY
LABCOLD

Robustly constructed from steel with a stainless steel interior this refrigerator from the **Labcold** Advanced pharmacy range is designed to provide close temperature control for temperature sensitive items in the pharmacy such as vaccines.

In addition, remote alarm contacts can be fitted as well as access ports. Contact us for bespoke inventory systems or specialised applications.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is hygienic stainless steel making it easy to clean and stock. The 8 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience are self closing.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Pharmacy Refrigerator is ideal for demanding environments.

Part Number	RPFG44043
Capacity (Litres)	1300
Style	Freestanding
Exterior Dimensions (HxWxD)	1985 x 1420 x 800mm
Energy Consumption	3.6 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	White polyester coated sheet steel
Shelves	8
Maximum shelf loading kg	70
Shelf dimensions (W x D)	500 x 620 mm
Door Hinging	Opens from centre outwards
Door Style	Glass
Internal Light	✓
Lockable Castors	✓
2 Door Locks	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK only)	✓



LABCOLD™

medical & scientific refrigeration



RPF05043 Solid Door Pharmacy & Vaccine Refrigerator*

Part Number	RPF05043
Capacity (Litres)	150
Style	Under Counter
Exterior Dimensions (HxWxD)	820 x 600 x 640 mm
Interior Dimensions (HxWxD)	660 x 545 x 355 mm
Energy Consumption	1.7 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	530 x 325 mm
Door Hinging	Right hand only
Door Style	Solid (Glass Door Model available RPF05043)
Internal Light	Glass Door Model only
Rear Castors	✓
Door Lock	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Auto defrost	✓
Fan Assisted Cooling	✓
Environmentally friendly insulation	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK Only)	✓
DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant	✓

Robustly constructed from Stainless Steel this refrigerator from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products.



The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.

Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.

* Glass door model available. Please order RPF05043.

01256 705 570

www.labcold.com



RPFR21043 Solid Door Advanced Pharmacy Refrigerator

CLAIM
YOUR
FREE
2
YEAR
WARRANTY
LABCOLD



Robustly constructed from steel with a stainless steel interior this refrigerator from the **Labcold** Advanced Pharmacy Range is designed to provide close temperature control for temperature sensitive medicines in the pharmacy.

In addition, remote alarm contacts can be fitted as well as access ports. Contact us for bespoke racking inventory systems or specialised applications.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is hygienic stainless steel making it easy to clean and stock. The 4 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience are self closing.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Pharmacy Refrigerator is ideal for demanding environments.

Part Number	RPFR21043
Capacity (Litres)	650
Style	Freestanding
Exterior Dimensions (HxWxD)	1985 x 710 x 800mm
Energy Consumption	2.43kWh/24h
Operating Temperature	+5°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	White polyester coated sheet steel
Shelves	4
Maximum shelf loading kg	70
Shelf dimensions (W x D)	500 x 620 mm
Door Hinging	Right side hinging only
Door Style	Solid
Internal Light	✓
Lockable Castors	✓
Door Lock	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK Only)	✓



RPFG21043 Glass Door Advanced Pharmacy Refrigerator

Part Number	RPFG21043
Capacity (Litres)	650
Style	Freestanding
Exterior Dimensions (HxWxD)	1985 x 710 x 800mm
Energy Consumption	2.43 kWh/24h
Operating Temperature	+5°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	White polyester coated sheet steel
Shelves	4
Maximum shelf loading kg	70
Shelf dimensions (W x D)	500 x 620 mm
Door Hinging	Right side hinging only
Door Style	Glass
Internal Light	✓
Lockable Castors	✓
Door Lock	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK only)	✓

Robustly constructed from steel with a stainless steel interior this refrigerator from the **Labcold** Advanced Pharmacy Range is designed to provide close temperature control for temperature sensitive medicines in the pharmacy.

In addition, remote alarm contacts can be fitted as well as access ports. Contact us for bespoke racking inventory systems or specialised applications.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is hygienic stainless steel making it easy to clean and stock. The 4 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience are self closing.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Pharmacy Refrigerator is ideal for demanding environments.



LABCOLD™

medical & scientific refrigeration



RPFR44043 Solid Door Advanced Pharmacy Refrigerator

Part Number	RPFR44043
Capacity (Litres)	1300
Style	Freestanding
Exterior Dimensions (HxWxD)	1985 x 1420 x 800mm
Energy Consumption	3.6 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	White polyester coated sheet steel
Shelves	8
Maximum shelf loading kg	70
Shelf dimensions (W x D)	500 x 620 mm
Door Hinging	Opens from centre outwards
Door Style	Solid
Internal Light	✓
Lockable Castors	✓
2 Door Locks	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Low GWP Refrigerant	✓
2 Year Warranty (UK Only)	✓

Robustly constructed from steel with a stainless steel interior this refrigerator from the **Labcold** Advanced pharmacy range is designed to provide close temperature control for temperature sensitive items in the pharmacy.

In addition, remote alarm contacts can be fitted as well as access ports. Contact us for bespoke racking inventory systems or specialised applications.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is hygienic stainless steel making it easy to clean and stock. The 8 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience are self closing.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Pharmacy Refrigerator is ideal for demanding environments.

01256 705 570

www.labcold.com

