



- 40°C Labcold Laboratory Freezers





Warning: Do not store open containers of volatile substances in this freezer.

Compact chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

RLHE0845 -40°C Sparkfree Chest Freezer

Part Number	RLHE0845
Capacity (Litres)	237
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1050 x 730mm
Interior Dimensions (HxWxD)	640 x 845 x 445mm
Weight	120kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Door Lock	✓
Feet	4
Energy Consumption	0.61 kWh/24h
Eco Friendly Refrigerant	✓



Warning: Do not store open containers of volatile substances in this freezer.

Compact chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

RLHE1145 -40°C Sparkfree Chest Freezer

Part Number	RLHE1145
Capacity (Litres)	314
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1300 x 730mm
Interior Dimensions (HxWxD)	665 x 1100 x 450mm
Weight	58kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
Door Open Alarm	✗
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Feet	4
Energy Consumption	0.71 kWh/24h
Eco Friendly Refrigerant	✓

RLHE1345 -40°C Sparkfree Chest Freezer



Warning: Do not store open containers of volatile substances in this freezer.

Robust chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

Part Number	RLHE1345
Capacity (Litres)	375
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1500 x 730mm
Interior Dimensions (HxWxD)	640 x 1295 x 445mm
Weight	68kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
Door Open Alarm	X
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Feet	4
Energy Consumption	0.8 kWh/24h
Eco Friendly Refrigerant	✓



RAVF1840 -40°C Laboratory Freezer

Part Number	RAVF1840
Capacity (Litres)	505
Style	Freestanding
Exterior Dimensions (HxWxD)	1702 x 770 x 725mm
Interior Dimensions (HxWxD)	1490 x 645 x 545mm
Energy Consumption	4.1 kWh/24h
Weight	86kg
Operating Temperature	-40°C
Interior Construction	White Enamel Steel
Exterior Construction	White Enamel Steel
Fixed Shelves	6
Door Type	Solid
Door Hinging	Right
Reverse Hinged	Optional
Defrost	Manual
Access Port	Optional
Open Door Alarm	X
High Temperature Alarm	✓
Low Temperature Alarm	✓
Remote Alarm Contacts	Optional
Digital Temperature Display	✓
Door Lock	✓
Feet	2 + 2 castors
Eco Friendly Refrigerant	✓

This **Labcold** freezer is constructed from enamelled mild steel with an enamel steel interior so is built to last. Ideal for any application which requires the storage of a large quantity of product within tight temperature tolerances. This freezer is manual defrost, putting you completely in control of the temperature of your samples.

**Looking for
temperature
mapping or
calibration?
Call the experts on
01256 705 570.**

Made to order at our UK factory, this freezer is specifically designed to withstand the rigours of professional use, and is backed by a comprehensive one year parts and labour warranty.

Fitted with a digital temperature display, you can be sure with a quick glance that the freezer is operating at the correct temperature. For extra peace of mind this freezer also has an audio/visual minimum and maximum temperature alarm and can have remote alarm contacts fitted at time of manufacture. For extra security this freezer comes complete with a door lock and two keys.