



- 80°C Labcold Upright Laboratory Freezers



Labcold -80°C Ultra Low Temperature Freezer LULT3585



Using tried and tested cascade refrigeration this quiet and compact Ultra Low Temperature freezer from Labcold features VIP insulation to improve energy usage and maximise internal space without increasing external size.

The controller displays the temperature and features high and low temperature and a door alarm plus it has battery back up so will function in the event of a mains failure.



Intuitive Control Panel

This freezer also has a door lock plus has an access port and remote alarm contacts and can fit up to 5 stainless steel drawers (optional extra) to ensure the contents can be organised.



For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.

Part Number	LULT3585
Capacity (Litres)	35
Style	Desk Top
Exterior Dimensions (HxWxD)	540 x 580 x 765 mm
Interior Dimensions (HxWxD)	280 x 425 x 300 mm
Weight	81kg
Operating Temperature	-80°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
VIP insulation	✓
Door Style	Solid
Door Hinging	Right
Open Door alarm	✓
High Temperature Alarm	✓
Mains Failure Alarm	✓
Digital Temperature Display	✓
Access Port	✓
Remote Alarm Contacts	✓
Key Operated Door Lock	✓
Battery Back up (Display & Alarms only)	✓
Drawers	5 (optional extra)
Inner Door	✗
Cascade Refrigeration	✓
FREE 5 year parts and 3 year labour warranty (UK only)	✓
Energy Consumption (kWh/24h)	4
Low GWP Refrigerant	✓



LABCOLD™

medical & scientific refrigeration



Using tried and tested cascade refrigeration this new under bench ULT freezer from Labcold also features VIP insulation to keep the footprint even smaller.



Intuitive Control Panel

For sample security the LULT80100 has a door lock plus features all the relevant alarms. It also comes complete with factory fitted remote alarm contacts and an access port.

For extra peace of mind, this freezer has a FREE 5 year parts warranty and a 3 year labour warranty.



**Need racking?
Why not ask for a
quote?**

Labcold -80°C Ultra Low Temperature Freezer LULT80100

Part Number	LULT80100
Capacity (Litres)	100
Style	Under bench
Exterior Dimensions (HxWxD)	800 x 950 x 735 mm
Interior Dimensions (HxWxD)	500 x 450 x 450 mm
Weight	120 kg
Operating Temperature	-80°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
VIP insulation	✓
Door Style	Solid
Door Hinging	Right
Open Door alarm	✓
High Temperature Alarm	✓
Mains Failure Alarm	✓
Digital Temperature Display	✓
Access Port	✓
Remote Alarm Contacts	✓
Key Operated Door Lock	✓
Battery Back up (Display & Alarms only)	✓
Shelves	1
Inner Door	✗
Rollers	✓
Cascade Refrigeration	✓
FREE 5 year parts and 3 year labour warranty (UK only)	✓
Energy Consumption (kWh/24h)	8.7
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



Labcold -50 to -86°C Ultra Low Temperature 390 litre Freezer

Part Number	LULTO400U
Capacity (Litres)	390
Style	Freestanding
Exterior Dimensions (HxWxD)	1990 x 840 x 750mm
Interior Dimensions (HxWxD)	1310 x 610 x 500mm
Weight	280kg
Operating Temperature	-50°C to -86°C
Interior Construction	Stainless Steel
Exterior Construction	Epoxy Painted Steel
Door Style	Solid
Door Hinging	Right
Remote Alarm contacts	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Mains Failure Alarm	✓
Door open alarm	✓
Low Battery Alarm	✓
Condenser High Temperature alarm	✓
Digital Temperature Display	✓
Key Operated Door Lock	✗
Password Protected Temperature Settings	✓
Battery Back up (Display and Alarms)	✓
Shelves	3
Inner Doors	4
Castors	4
Access Port	✓
Eco Friendly Refrigerant	✓
5 year warranty 3 x years parts and labour 2 x years parts only	✓

Labcold Ultra Low Temperature freezers use the tried and trusted cascade refrigeration and have been designed to occupy less than 1m² of floor space.



Complete with an easy to use digital controller conveniently situated on the front of the freezer, the temperature settings can be password protected and it always displays the set point for extra peace of mind. It also features alarm battery

back -up so should there be a mains power failure the freezer will alarm immediately to give you as much time as possible to deal with the emergency.

Backed by a 5 year parts and 3 year labour warranty, this Ultra Low Temperature freezer from Labcold is designed to be a robust and reliable place to store your samples.

We can also supply racking, why not ask for details at time of order?



01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



Labcold -50 to -86°C Ultra Low Temperature 500 litre Freezer

Part Number	LULT0500U
Capacity (Litres)	500
Style	Freestanding
Exterior Dimensions (HxWxD)	1990 x 840 x 880mm
Interior Dimensions (HxWxD)	1310 x 600 x 630mm
Weight	320kg
Operating Temperature	-50°C to -86°C
Interior Construction	Stainless Steel
Exterior Construction	Epoxy Painted Steel
Door Style	Solid
Door Hinging	Right
Remote Alarm contacts	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Mains Failure Alarm	✓
Door Open Alarm	✓
Low Battery Alarm	✓
Condenser High Temperature alarm	✓
Digital Temperature Display	✓
Key Operated Door Lock	✗
Password Protected Temperature Settings	✓
Battery Back up (Display and Alarms)	✓
Shelves	3
Inner Doors	4
Castors	4
Access Port	✓
Eco Friendly Refrigerant	✓
5 year warranty 3 x years parts and labour 2 x years parts only	✓

Labcold Ultra Low Temperature freezers use the tried and trusted cascade refrigeration and have been designed to occupy less than 1m² of floor space.



Complete with an easy to use digital controller conveniently situated on the front of the freezer, the temperature settings can be password protected and it always displays the set point for extra peace of mind. It also features alarm battery

back -up so should there be a mains power failure the freezer will alarm immediately to give you as much time as possible to deal with the emergency.

Backed by a 5 year parts and 3 year labour warranty, this Ultra Low Temperature freezer from Labcold is designed to be a robust and reliable place to store your samples.

We can also supply racking, why not ask for details at time of order?



01256 705 570

www.labcold.com

