



Sparkfree Freezers





LABCOLD™

medical & scientific refrigeration



RLVF0217 Bench Top Sparkfree Freezer

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



CLAIM
YOUR
FREE
YEAR
WARRANTY

For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.

Part Number	RLVF0217
Capacity (Litres)	66
Style	Bench Top
Exterior Dimensions (HxWxD)	740 x 450 x 510mm
Interior Dimensions (HxWxD)	590 x 333 x 300mm
Weight	34kg
Energy Consumption	0.19kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Shelves	3
Shelf Dimensions (w x d)	2 of 309 x 275mm 1 of 309 x 150mm
Door Lock	✓
Door Open Alarm	✓
Feet	4
Access Port	✓
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



RLVF0417 Sparkfree Freezer

Part Number	RLVF0417
Capacity (Litres)	124
Style	Under Bench
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Interior Dimensions (HxWxD)	680 x 460 x 410mm
Weight	45kg
Energy Consumption	0.8kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Powder Coated Steel
Exterior Construction	Powder Coated Steel
Door Style	Self Closing
Door Hinging	Right (can be reversed please ask at time of order)
Shelves	3
Shelf Dimensions (w x d)	2 of 460 x 407mm 1 of 460 x 240mm
Door Lock	✓
Door Open Alarm	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Access Port	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



RLVF1517 Freestanding Sparkfree Freezer

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm.

Part Number	RLVF1517
Capacity (Litres)	406
Style	Free standing
Exterior Dimensions (HxWxD)	1875 x 600 x 700mm
Interior Dimensions (HxWxD)	1690 x 440 x 500mm
Net Weight	97kg
Energy Consumption	3.16kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Painted Aluminium
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	413mm x 438mm
Fan Circulated	✓
Rollers	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Access Port	Optional
Automatic Defrost	✗
2 Year Warranty (UK only)	✓

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



RLCF0720 Sparkfree Chest Freezer

Part Number	RLCF0720
Capacity (Litres)	215
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 725 x 655mm (including Handle)
Interior Dimensions (HxWxD)	697 x 606 x 536mm
Weight	46kg
Energy Consumption	0.71 kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	1
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
High Temperature Alarm	Visual
Defrost drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with a storage basket to assist with organisation and storage. Please ask at time of order for additional baskets. This chest freezer has manual defrost for greater operator control.

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration



RLCF1520 Sparkfree Chest Freezer

Part Number	RLCF1520
Capacity (Litres)	527
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1505 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	697 x 1386 x 536mm
Weight	67kg
Energy Consumption	1.3kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Lock	✓
Baskets	3
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security. This model comes with 3 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets or racking.

This chest freezer has manual defrost for greater operator control.

01256 705 570

www.labcold.com





LABCOLD™

medical & scientific refrigeration

RLCF2120 Sparkfree Chest Freezer



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this freezer.

Did you know that 25% of chemical accidents are caused by improper storage?

Source: US Chemical Safety and Laboratory Investigation Board

Designed specifically for laboratory safety with a sparkfree interior and factory fitted with a door lock for optimum security.

This model comes with 3 storage baskets to assist with organisation and storage. Please ask at time of order for additional baskets.

This chest freezer has manual defrost for greater operator control and a defrost drain to assist with getting rid of any defrost liquids quickly plus it comes with a 1 year UK parts and labour warranty.

Part Number	RLCF2120
Capacity (Litres)	607
Style	Chest Freezer
Exterior Dimensions (HxWxD)	865 x 1705 x 655mm (Including Handle)
Interior Dimensions (HxWxD)	697 x 1586 x 536mm
Weight	77kg
Energy Consumption	1.52kWh/24hr
Operating Temperature	-18°C to -25°C
Factory Set Point	-20°C
Interior Construction	Aluminium
Exterior Construction	White Enamel
Lid Style	Solid
Door Hinging	Back (viewed from front)
Door Lock	✓
Baskets	3
Basket Dimension (HxWxD)	200 x 250 x 500mm
Feet	4 (adjustable)
Mains Indicator Light	✓
Defrost Drain	✓
Manual Defrost	✓
Eco Friendly Refrigerant	✓

01256 705 570

www.labcold.com

